

E.Z.N.A.® DNA/RNA Clean-Up System Selection Guide

	E.Z.N.A.® Clean-Up Product	Size Range	Max Recovery	Elution volume	Format	Processing
PCR Clean-Up	Cycle-Pure Kit	100bp-10kb	85%	30-50µL	Spin Columns	spin/vacuum
	E-Z 96® Cycle-Pure Kit	100bp-10kb	85%	80-100µL	96-well plates	spin/vacuum
	MicroElute® Cycle-Pure Kit	100bp-10kb	85%	10-15µL	Spin Columns	spin/vacuum
	E-Z 96® MagBind® Cycle-Pure Kit	100bp-10kb	90-95%	20-50µL	96-well plates	magnet
Gel Extraction	Gel Extraction Kit	70bp-20kb	85%	30-50µL	Spin Columns	spin/vacuum
	MicroElute® Gel Extraction Kit	70bp-20kb	85%	10-15µL	Spin Columns	spin/vacuum
	Ultra-Sep® Gel Extraction Kit	60bp-200bp 200bp-1000bp 1-5kb 25-50kb	< 60% 85% 85% 80%	≥20µL	Silica particles	spin
	Poly-Gel DNA Extraction Kit	100bp-10kb	n/a	30-50µL	Spin Columns	spin
	Mag-Bind® Poly-Gel DNA Extraction Kit	≤100bp	n/a	20-50µL	paramagnetic beads	magnet
DNA Clean-Up	MicroElute® DNA Clean-Up Kit	100bp-10kb	85%	10-15µL	Spin Columns	spin/vacuum
	DNA Probe Purification Kit	100bp-10kb	85%	30-50µL	Spin Columns	spin
Nucleotide Removal	Mag-Bind® Oligonucleotide Purification Kit	10bp-100bp	60-75%	20-50µL	paramagnetic beads	magnet
	E-Z 96® Mag-Bind® Oligonucleotide Purification Kit	10bp-100bp	60-75%	30-50µL	paramagnetic beads	magnet
Dye-Terminator Removal	Ultra-Sep® Dye Terminator Removal Kit	n/a	n/a	n/a	prehydrated columns	spin
	E-Z 96® Ultra-Sep® Dye Terminator Removal Kit	n/a	n/a	n/a	prehydrated columns	spin
	E-Z 96® Mag-Bind® Dye Terminator Removal Kit	n/a	n/a	n/a	prehydrated 96-Well plate	magnet
RNA Cleanup	Poly-Gel RNA Extraction Kit	n/a	n/a	30-50µL	spin columns	spin
	RNA Probe Purification Kit	n/a	85-90%	30-50µL	spin columns	spin
	MicroElute® RNA Clean-Up Kit	n/a	85-90%	10-15µL	spin columns	spin
	Mag-Bind® RNA Clean-Up Kit	n/a	85-90%	5-10µL	paramagnetic	magnet